



VITALCOR INC.

MEDICAL TECHNOLOGY IN MOTION

Gemini Plus Surgical Headlight
Product Number 847211




The Gemini Plus Headlight. The spot size has been increased to 140 mm and the improved optics creates a brighter light regardless of the light source. The encapsulated ultra-low-profile design does not interfere with surgical loupes and allows surgeons to work closely without bumping headlights. The headlight weighs 8 ounces and features two easy –to-use controls for spot size and lens adjustment. The sweatbands are removable and replaceable.

Available with Camera Soon



Spot adjustment knob



Lens adjustment knob



Strap keeps cable in place and secure



Size adjustment knobs on the top and back of headlight



Works well with Loupes



Eliminates the "ramming" Effect



Gemini Plus Headlight
Product No. 847211
Fiber Optic Cable Product No. 847310-Clear
Cables available in different colors