



Coronary Artery Perfusion Cannulae with Balloon

Coronary Perfusion/Cardioplegia Infusion Cannulae

Self Inflating PVC Balloon (No Latex)

Female Luer Lock Fitting

Balloon Protection With Unique Packaging. Longer Shelf Life Of Three Years



Coronary Artery Perfusion Cannulae with Balloon

Total Length 11.4 in./29 cm

Balloon Size	Distal Cannulae Tip Size	Catalog number	
		Straight	Right Angle
4.0 mm	2.1 mm	315804	415804
5.0 mm	2.1 mm	315805*	415805*
6.0 mm	3.0 mm	315806	415806
7.0 mm	3.0 mm	315807	415807
8.0 mm	3.5 mm	315808	415808

* New Size Balloon
5 Cannulae per Box. Ethylene Oxide Sterilized in Individual Pouch.
Single Use Only. United States Origin & Manufacture.

VITALCOR CE
MEDICAL TECHNOLOGY IN MOTION

100 East Chestnut Avenue • Westmont, Illinois 60559
Phone: In Illinois (630) 325-5500 • Outside Illinois (800) 874-8358
Fax: (630) 325-0257 • <http://www.vitalcor.com>

Consult product insert sheet accompanying individual products for instruction on use, indications, contraindications, warnings and precautions
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician